Eja:<on





Category


Topical gel for erectile dysfunction treatment

How to Use

Apply the contents of one tube to the head of the penis

When to Use

Just before sex



Indication



A clinically proven product for the treatment of erectile dysfunction in men aged 18 years and over

Mechanism of Action



A unique physical cooling then warming effect following the application to the head of the penis and massaging for 15 seconds



Stimulation of penile nerve endings, increased blood flow and getting an erection within 10 minutes







Ingredients	
Ethanol (35%)	Propylene glycol
Glycerin	Carbomer
Potassium hydroxide	Aqua
Triethanolamine	

Packaging

4 single-dose tubes per box



References

- 1- https://www.drugs.com/medical-answers/what-eroxon-3574126/
- 2- https://www.medscape.com/viewarticle/otc-solution-erectile-dysfunction 2024a10008ms.
- 3- Ejaxon product leaflet, Espad Pharmed company.